

## INTRODUCTION




Polymeric microfilm matrix

- ❖ Biofilms are bacteria encased in extracellular polymeric substance (EPS)
- ❖ Biofilms exist in 60% of chronic wounds, 6% of acute wounds<sup>1</sup>
- ❖ Biofilms are diffusive barriers that make bacteria 1000 times more resistant to antibiotics/antimicrobials than planktonic bacteria<sup>2,3</sup>

**Goal:** To identify combination of **gallium (antibiofilm agent)** and **silver (antimicrobial)** that works synergistically within a **polymeric microfilm to disperse biofilms and kill bacteria.**

## INNOVATION


**Strategy:** Wound surface modification with **synthetic polymer microfilm that releases silver and gallium**




Microfilm consists of two layers

- ❖ Water soluble polymer cast with or without gallium nitrate
- ❖ Polyelectrolyte multilayers (PEMs) matrix containing ionic silver and gallium and metallic silver

## Antibiofilm Gallium and Antimicrobial Silver




- ❖ Gallium **replaces iron needed for metabolism** due to its chemical similarities<sup>4,5</sup>
- ❖ Gallium cannot be reduced under physiological conditions<sup>4,5</sup>
- ❖ Gallium **blocks redox-driven biological processes of metabolism** that reduces bacterial growth and biofilm formation<sup>4,5</sup>
- ❖ 10 µg/cm<sup>2</sup> silver in PEMs
- ❖ 5-105 µg/cm<sup>2</sup> gallium in microfilm



Films extracted in simulated wound fluid (SWF, 1x MEM with 10% FBS) at 37°C

- ❖ Burst release of gallium as polymer cast dissolves
- ❖ Sustained release of silver as silver nanoparticles dissolve




MTT Cytotoxicity Test

- ❖ Gallium in microfilm extracted in SWF at 37°C for 24 hours
- ❖ Extract dosed on L929 mouse fibroblasts for 24 hours at 37°C
- ❖ 62% cell viability at 70 µg/cm<sup>2</sup> gallium
- ❖ Dose limits of gallium 0-125 µg/cm<sup>2</sup>

## METHODS

Samples evaluated in-vitro against 48-hr old mature biofilms grown on cotton gauze




Post treatment, biofilms supported on gauze were rinsed free of planktonic bacteria and transferred to clean well plates.

- ❖ Samples stained with crystal violet and absorbance (@395 nm wavelength) read as a representation of biofilm mass
- ❖ Samples homogenized and serial plated to quantify *P. aeruginosa* colony forming units (CFU)


## RESULTS

❖ With no silver in microfilm




- ❖ Gallium alone is ineffective at reducing bacteria contained within biofilm
  - Limited to **0.6 ± 0.4 log<sub>10</sub> reduction** of *P. aeruginosa* with **70 µg/cm<sup>2</sup> gallium in microfilm**
- ❖ Gallium alone is ineffective at dispersing biofilm mass
  - Achieved **less than 25% biofilm mass dispersal** with up to **125 µg/cm<sup>2</sup> gallium in microfilm**

❖ With 75 µg/cm<sup>2</sup> gallium in microfilm, varied silver quantity within PEMs layer



✓ With **addition of silver** to 75 µg/cm<sup>2</sup> of gallium, **>4.5 ± 0.6 log<sub>10</sub> CFU reduction to *P. aeruginosa* can be achieved**

❖ With 10 µg/cm<sup>2</sup> of silver in PEMs, varied gallium level in microfilm.



✓ 10 µg/cm<sup>2</sup> silver and 5 µg/cm<sup>2</sup> gallium in microfilm achieves **>80% biofilm mass dispersal**

✓ When microfilm contains more than 35 µg/cm<sup>2</sup> gallium the prototype **reduces more viable bacteria encased in biofilms than commercial silver products**

## CONCLUSIONS

Gallium has limited effectiveness against bacteria contained in biofilms at nontoxic doses

Silver has limited effectiveness against bacteria contained in biofilms at nontoxic doses

Silver and gallium together can achieve >4log<sub>10</sub> reduction and disperse 80% biofilm mass at benign doses

## REFERENCES

- (1) James GA, et al. Biofilms in chronic wounds. Wound repair and regeneration : Official Publication of the Wound Healing Society [and] the European Tissue Repair Society. 2008;16(1):37-44.
- (2) Rhoads DD, et al. Percival SL. Biofilms in wounds: management strategies. *Journal of Wound Care*. 2008;17(11):502-8.
- (3) Stewart PS, et al. Antibiotic resistance of bacteria in biofilms. *The Lancet*. 2001;358(9276):135-8.
- (4) Herron, Maggie, et al. "Gallium-Loaded Dissolvable Microfilm Constructs That Provide Sustained Release of Ga<sup>3+</sup> for Management of Biofilms." *Advanced Healthcare Materials*, vol. 4, no. 18, 2015, pp. 2849–2859., doi:10.1002/adhm.201500599.
- (5) Herron, Maggie, et al. "Interfacial Stacks of Polymeric Nanofilms on Soft Biological Surfaces That Release Multiple Agents." *ACS Applied Materials & Interfaces*, 8(40) 2016:26541.