

# Effective Antibiofilm Bioresorbable Microfilm™ Matrix Containing Gallium and Silver

Gaurav Pranami, PhD<sup>1</sup>, Eric C. Crawford<sup>1</sup>, Ankit Agarwal, PhD<sup>1</sup>, Nicholas L. Abbott, PhD<sup>2</sup>, Patricia R. Kierski<sup>3</sup>, Charles J. Czuprynski, PhD<sup>3</sup> and Jonathan F. McAnulty, DVM, PhD<sup>3</sup>, Michael J. Schurr, MD<sup>4</sup>

<sup>1</sup>Imbed Biosciences Inc., Madison, WI, <sup>2</sup>Department of Chemical and Biological Engineering, University of Wisconsin-Madison, <sup>3</sup>Department of Surgical Sciences, School of Veterinary Medicine, University of Wisconsin-Madison, <sup>4</sup>Department of Surgery, Mission Hospitals, Asheville, NC

## INTRODUCTION

The health care costs associated with treatment of chronic and burn wounds exceeds \$25 billion annually in the U.S., involving 8.6 million patients (1). Biofilms play a significant role in rendering a wound chronic as evidenced by James *et al.* who showed that **60% of chronic wound debridement samples contained biofilm, whereas only 6% of acute wound biopsies contained biofilm (2).**

Biofilms are structured communities of bacteria encased and protected in extracellular polymer matrix (EPM), which comprises 80% volume (3).

Due to restricted penetration of antimicrobial agents through EPM, presence of persister cells and other mechanisms, **bacteria in biofilms are 1000 times more difficult to kill than planktonic bacteria (4,5).**

Many wounds have complex surfaces and debridement can be challenging, leaving biofilm fragments that remain resistant to antimicrobial therapy and act as a nidus for recrudescence of biofilms. **There is currently no commercially available wound dressing that is clinically indicated for dispersal of biofilms in a wound bed.** Commercial antimicrobial dressings have been shown to only be marginally effective against biofilms delivering 1-2 Log<sub>10</sub> CFU reduction (6).

## INNOVATION

Several studies have shown that gallium (Ga) prevents bacterial growth, biofilm formation and disperses established biofilms (7,8). However, at non-toxic concentrations Ga<sup>3+</sup> is not effective at killing bacteria. On the other hand, Ag<sup>+</sup> is a broad spectrum antimicrobial effective in killing planktonic bacteria but not the bacteria encased in biofilms.

**We are the first to demonstrate that synergy between non-cytotoxic levels of Ag<sup>+</sup> and Ga<sup>3+</sup> disperses biofilm and kills the bacteria encased in it (9).** Specifically, we determined that non-cytotoxic concentrations of 1.5 µg/ml Ag<sup>+</sup> (14 µM) and 10 µg/ml Ga<sup>3+</sup> (143 µM), when used independently, lowered biofilm CFU by only 20% and 25%, respectively. However, an aqueous solution containing 1.5 µg/ml of Ag<sup>+</sup> and 10 µg/ml of Ga<sup>3+</sup> lowered biofilm CFU by 85%.

## REFERENCES

- Chandan *et al.* Human skin wounds: A major and snowballing threat to public health and the economy. *Wound Repair and Regeneration*. 2009;17(6):763-71.
- James *et al.* Biofilms in chronic wounds. *Wound repair and regeneration: official publication of the Wound Healing Society [and] the European Tissue Repair Society*. 2008;16(1):37-44.
- Metcalf *et al.* Biofilm delays wound healing: A review of the evidence. *Burns & Trauma*. 2013; 1(1):5-12.
- del Pozo *et al.* The challenge of treating biofilm-associated bacterial infections. *Clinical Pharmacology & Therapeutics*. 2007; 82(2):204-9.
- Rhoads *et al.* Biofilms in wounds: management strategies. *Journal of Wound Care*. 2008;17(11):502-8.
- Kostenko *et al.* Impact of silver-containing wound dressings on bacterial biofilm viability and susceptibility to antibiotics during prolonged treatment. *Antimicrobial agents and chemotherapy*. 2010;54(12):5120-31.
- Kaneko *et al.* The transition metal gallium disrupts *Pseudomonas aeruginosa* iron metabolism and has antimicrobial and antibiofilm activity. *The Journal of Clinical Investigation*. 2007;117(4):877-88.
- Zhu *et al.* Pre-treatment with EDTA-gallium prevents the formation of biofilms on surfaces. *Experimental and therapeutic medicine*. 2013;5(4):1001-4.
- Herron *et al.* Interfacial Stacks of Polymeric Nanofilms on Soft Biological Surfaces that Release Multiple Agents. *ACS Applied Materials & Interfaces*. 2016;8(40):26541-51.

## ACKNOWLEDGEMENT

This research was supported by an NIH SBIR Phase II grant from NIAMS (5R43AR061913-02A1). AA, MJS, CJC and JFM possess financial interests in *Wound Engineering LLC* and/or *Imbed Biosciences Inc.*, for-profit organizations that have filed patent applications and/or are commercializing aspects of the work reported in this presentation.

## METHODS

Microfilm matrix was prepared as follows,

An interpenetrated 300 nm thick nanofilm was assembled by alternating deposition of oppositely charged polyelectrolytes (polyacid and polyamine). Post-assembly, nanofilm was impregnated with Ag<sup>+</sup> ions via ion exchange and Ag nanoparticles (5-20 nm) were formed in-situ via reduction. Nanofilm was further impregnated with Ga<sup>3+</sup> ions via ion exchange, as needed. A bioresorbable 20 µm thick polyvinyl alcohol hydrogel w/wo Ga<sup>3+</sup> ions was laminated to nanofilm.




Fig. (A) Schematic of bioresorbable matrix containing silver and gallium, (B) TEM image showing silver nanoparticles embedded into nanofilm, (C) Bioresorbable matrix

## P. aeruginosa biofilm growth,

Robust biofilms of *P. aeruginosa* were grown on 15 mm squares of Kendall Curity Gauze Bandage (Tyco Healthcare, KDL1143) by plating 800 µl suspension of *P. aeruginosa* (10<sup>8</sup> CFU/ml) for 48 h at 30 °C. Post incubation, gauze specimen was rinsed 3x with sterile 1x PBS to remove loosely adhering planktonic bacteria. The presence of biofilm was confirmed with crystal violet (CV) staining which binds specifically to EPM. Biofilm mass was quantified as absorbance of CV at 595 nm using a spectrophotometer.




Fig. CV absorbance of 48 h *P. aeruginosa* biofilm relative to background confirms robust biofilm growth on gauze

## Antibiofilm testing,

Gauze specimen supporting biofilm were placed on tryptic soy agar contact plate (w 5% sheep blood), hydrated with 175 µl of 1x PBS and treated with 20 mm x 20 mm test matrix for 24 h at in a humidified incubator at 30 °C.




Fig. (A) Untreated gauze supporting biofilm (control), (B) Microfilm matrix applied over gauze supporting biofilm




Fig. Schematic showing that unlike (A) conventional silver dressings, (B) Microfilm matrix adheres intimately to wound bed resulting in a targeted and localized action


Post treatment, gauze specimens were transferred to 24 well plates and rinsed 3x with sterile PBS. One set of specimens was homogenized and serial plated to determine biofilm CFU, and another set was treated with CV to determine total biofilm mass.

## RESULTS


Commercial antimicrobial dressings Silverlon<sup>®</sup> and Aquacel<sup>®</sup> Ag Extra containing high level of silver, 5460 and 120 µg/cm<sup>2</sup> respectively, are not only ineffective against *P. aeruginosa* biofilms *in-vitro* but also cytotoxic

Silverlon and Aquacel Ag Extra showed no reduction in *P. aeruginosa* biofilm CFU over 24 h relative to untreated control.


Dressings were extracted in growth media for 24 h at 37°C; extracts and their dilutions were incubated with L-929 mouse fibroblasts cells in growth media (5% FBS) for 24 h; Cell viability assessed by NAMS (Northwood, OH) using MTT cell metabolic activity assay. At least 70% cell viability indicates non-cytotoxicity.



Microfilm matrix containing only 10 µg/cm<sup>2</sup> of silver kills 4 Log<sub>10</sub> CFU of planktonic bacteria per ISO 22916 for up to 3 days without causing cytotoxicity



Microfilm matrix pairing 10 µg/cm<sup>2</sup> of silver with 70-165 µg/cm<sup>2</sup> of gallium synergistically kills 4-9 Log<sub>10</sub> CFU of *P. aeruginosa* biofilm bacteria and disperses > 90% biofilm mass



## CONCLUSION

Microfilm matrix amplifies the efficacy of silver and gallium ions 1000 times compared to topical aqueous formulations for combating biofilms.

## FUTURE WORK

We will evaluate the efficacy of microfilm matrix containing Ag and Ga in the healing of splinted mice wounds inoculated with mixed species (*S. aureus* and *P. aeruginosa*) biofilm, which impairs wound healing (Fig. on right), as shown by Dr. Jonathan F. McAnulty.




Fig. Healing of wounds covered for 3 days with gauze with a *S. aureus* biofilm. Wounds were rinsed on day 3 and then onwards treated with fresh gauze. Biofilm wounds had significantly slower healing over 21 days compared to control wounds with no biofilms (\* p<0.05; \*\* p<0.0001).



Fig. (Top) A splinted full-thickness wound-model in mice. (Bottom) Microfilm Matrix placed on a contaminated moist wound.